

# A new generation of environmental sampling devices for the food and pharmaceutical industries.



Coliform Bacteria



Salmonella

## NRS II Transwab®




## Reliable microbiological sampling of critical surfaces

### Neutralises trace disinfectant

NRS II Transwab® is MWE's name for its range of swab based environmental sampling devices for the food, pharmaceutical, biotechnology and cosmetic industries. There can also be applications in healthcare including the monitoring of infection control measures.

All NRS II Transwab® devices feature MWE's leak proof labelled self standing screw-cap tube made from shatterproof polypropylene, with a high visibility blue shaft swab attached to the cap. The swab features a rayon bud that can remain immersed in liquid, yet retains MWE's high standard for absorption, survival and release of microorganisms. The tubes are prefilled with the specified volume of solution and there is a choice of NRS, MWE's classic neutralising rinse solution, and a range of alternative buffers and media appropriate for particular applications.



## NRS (Neutralising Rinse Solution)

NRS Medium contains lecithin, polysorbate 80, and sodium thiosulphate in a peptone phosphate buffer. The formulation conforms to ISO 18593 and will neutralise most disinfectants used in the food industry, including those based on chlorine, peroxygen compounds, amphoterics, biguanides, and glutaraldehyde. This is a universal neutralising solution suitable for testing most disinfected areas within the food, cosmetic and pharmaceutical industries. Precise fill volumes allow accurate quantitative assessment of contamination levels.

## Buffered Peptone Water

Buffered Peptone Water contains peptone and sodium chloride, with a phosphate buffer, and is used for the recovery of Salmonella, especially in environments where they may have been sub-lethally injured.

## Butterfield Buffer

Buffer Butterfield Buffer (also called Butterfield's Solution or Butterfield's phosphate buffered dilution water) contains potassium phosphate as specified by APHA and FDA, and is used in standard methods for the enumeration of bacteria and fungi in foods, water and pharmaceutical products. The standard formulation has been modified by the inclusion of a low level of peptone to reduce osmotic shock, and polysorbate 80 as a surfactant to assist the sampling process.

## MRD (Maximum Recovery Diluent)

Also called Peptone Saline. Maximum Recovery Diluent is an isotonic and protective medium for maximum recovery of microorganisms from environmental and food sources. It is also the recommended diluent for the sampling of carcasses in the European Union.

## Tryptic Soy Broth

Tryptic Soy Broth is a general purpose recovery medium with excellent growth properties for aerobes, anaerobes, and some fungi. It is widely used for sterility testing, particularly in pharmaceutical manufacturing.

## Lethen Broth

Lethen media are used to determine the bactericidal efficacy of quaternary ammonium based disinfectants.



Lethen broth is a growth medium recommended by the FDA for use in the microbiological testing of cosmetics, and the formulation is as described in the FDA Bacteriological Analytical Manual. Polysorbate 80 and lecithin are included to partially neutralise the preservatives commonly found in cosmetics.

## Effective sampling of food contact surfaces after cleaning

### NRS without Sodium Thiosulphate

Some testing methods (such as Petrifilm®) require medium without thiosulphate.







### Choice of fill

NRS medium and variants are available in a choice of 3 fill sizes

- 10ml and 5ml are used with standard and filtration methods.
- 1ml can be used directly with Petrifilm®, and other direct pour plating techniques



#### Ordering info:

Part Number	Description	Quantity
 MW784	NRSII™ Transwab® 5ml fill, Blue	50
 MW785	NRSII™ Transwab® with 10ml fill, Blue	50
 MW786	NRSII™ Transwab® with 1ml fill, Blue	50
 MW792	Lethen Broth 10ml	50
 MW793	Lethen Broth 5ml	50
 MW794	Lethen Broth 1ml for Petrifilm as 3M Quick Swab	50

Distributed by



Free Call: 1800 723 405 - Ph: (03) 9845 0300 - [sales@pacificlab.com.au](mailto:sales@pacificlab.com.au)



[www.pacificlab.com.au](http://www.pacificlab.com.au) 