


Introducing PLP Cutting-Edge Cell Culture Plasticware High Quality, Transparency, and Value for Money



Highly Transparent: Our plasticware ensures clear visibility of your cell cultures, making it easy to observe growth

Premium Quality: Our cell culture plasticware is expertly crafted from Polystyrene (PS), providing optimal conditions for cell growth


Cost-Effective: Despite its high quality, our plasticware remains budget-friendly, offering exceptional value for your research needs

Versatile: Suitable for a wide range of applications

Great value

Cell Culture Dishes

- Optically clear Polystyrene for Microscopy
- Tissue Culture treated surface
- Unique lid for effective gas exchange
- Available in bulk or individually packaged
- RNase, DNase and Endotoxin Free
- Sterile



Cell Culture Dishes	Description	Packaging	Pack Size
Bulk Packed			
CCD-35	Cell Culture Dish, 35mm, TC treated, Sterile	20pc/Bag, 50 Bags/Carton	1000
CCD-60	Cell Culture Dish, 60mm, TC treated, Sterile	20pc/Bag, 25 Bags/Carton	500
CCD-100	Cell Culture Dish, 100mm, TC treated, Sterile	20pc/Bag, 25 Bags/Carton	500
CCD-150	Cell Culture Dish, 150mm, TC treated, Sterile	5pc/Bag, 6 Bags/Carton	30
Individually Packed			
CCD-35H	Cell Culture Dish, 35mm, TC treated, Sterile	Individually Packed	800
CCD-60H	Cell Culture Dish, 60mm, TC treated, Sterile	Individually Packed	200
CCD-100H	Cell Culture Dish, 100mm, TC treated, Sterile	Individually Packed	200
CCD-150H	Cell Culture Dish, 150mm, TC treated, Sterile	Individually Packed	60

Cell Culture Plates

- Designed to prevent cross contamination with effective gas exchange lid
- Flat bottom
- Tissue culture treated surface
- RNase, DNase and Endotoxin Free
- Sterile



Cell Culture Plates	Description	Packaging	Pack Size
CCP-6H	Cell Culture Plate, 6 Wells, Flat Bottom, TC treated with Lid, Sterile	Individually Packed	50
CCP-12H	Cell Culture Plate, 12 Wells, Flat Bottom, TC treated with Lid, Sterile	Individually Packed	50
CCP-24H	Cell Culture Plate, 24 Wells, Flat Bottom, TC treated with Lid, Sterile	Individually Packed	50
CCP-48H	Cell Culture Plate, 48 Wells, Flat Bottom, TC treated with Lid, Sterile	Individually Packed	50
CCP-96H	Cell Culture Plate, 96 Wells, Flat Bottom, TC treated with Lid, Sterile	Individually Packed	100

Cell Culture Flasks

- Vented Cap
- Tissue Culture treated surface
- Angled wide neck for easy access for serological pipettes and cell scrapers
- Stackable
- Available in bulk pack or individually packed
- RNase, DNase and Endotoxin Free
- Sterile



Cell Culture Flasks	Description	Packaging	Pack Size
Bulk Packed			
CCF-25F	Cell Culture Flask, 25cm ² Surface Area, Angled Neck, Vented Cap	10 Pc/Bag, 20 Bags/Carton	200
CCF-75F	Cell Culture Flask, 75cm ² Surface Area, Angled Neck, Vented Cap	5 Pc/Bag, 20 Bags/Carton	100
CCF-175F	Cell Culture Flask, 175cm ² Surface Area, Angled Neck, Vented Cap	5 Pc/Bag, 6 Bags/Carton	30
Individually Packed			
CCF-T25FH	Cell Culture Flask, 25cm ² Surface Area, Angled Neck, Vented Cap	Individually Packaged	100
CCF-T75FH	Cell Culture Flask, 75cm ² Surface Area, Angled Neck, Vented Cap	Individually Packaged	100
CCF-T175FH	Cell Culture Flask, 175cm ² Surface Area, Angled Neck, Vented Cap	Individually Packaged	60

Introductory prices are subject to change without notice. Prices are as quoted. All prices are in AUD and GST is applicable on all purchases. Local delivery fees apply. Further product information and specifications including procurement, Terms and Conditions of Sale can be found at www.pacificlab.com.au or by contacting sales@pacificlab.com.au. Images displayed on this publication are a graphic representation of the product only and may vary from the actual item itself.