



BenchTop Essentials SALE

Special offers on everyday lab equipment

Visit pacificlub.com.au

Please refer to the reverse page for terms and conditions.



Horizontal & Gyating Shaker with Remote Control

\$450^{.*}

ID: PLP-SYC-FP200



- Ideal to be used inside incubators, refrigerators, etc.
- Split-Type, a flat cable connects the controller to the unit
- Speed Range:
 - Pendulum Mode: 5-100rpm
 - Horizontal Mode: 5-200rpm
- Radius: 10 degrees
- Tray Size: 317 x 245mm
- Max. Load: 2kg

Tissue Homogeniser High-Throughput

Rapid tissue disruption using high-speed bead-based homogenisation. Designed for reliable, high-throughput processing of multiple samples.

- Room temperature 5°C - 40°C
- Relative humidity max. 80%
- High-efficiency sample processing
- Robust linear guide rail system
- Pre-cooling adapter module
- Flexible adapter options
- Compact, modern design

ID: PLP-SWE-C6

\$5,750^{.*}



Servicebio			
Sample	liver		
OP timer	20 s	OP speed	60 Hz
Pause timer	0 s	OP time	
		Pause time	
		OP times	




Hotplate Stirrer with Large Ceramic Plate

- Heating and stirring with intelligent speeding function
- Stirring volume capacity up to 10 litres
- Temperature sensor PT100 with 440mm rod and clamp
- LED display for temperature and speed range
- Automatic stirring bar homing, magnetic stir bar area
- Anti-scalding design, flame-retardant
- Over temperature protection

ID: PLP-SMS-200A

\$525^{.*}



Horizontal Shaker with Spring Clamp Platform



- Spring grid design
- Stable sample placement
- Automatic levelling
- LED display (speed/time)
- Speed Range: 5-200rpm
- Radius of Gyration: 15mm
- Max. Load: 5kg



\$550.*
.00

ID: PLP-SYC-WP200

Centrifuge Plate Spinner for 2 x 96 PCR plates



- Designed for PCR plates
- Emergency stop
- Brushless motor
- Large openings



ID: PLP-SMP-2

\$575.*
.00

Compact Laboratory Vacuum Aspirator



- Two liquid control methods: Foot pedal and knob switch
- Adjustable speed suction
- Max. Suction speed: 1.6L/min
- Air filtration: 0.45 µm
- Hydrophobic filtration membrane
- Pipetting bottle volume: 2L



\$550.*
.00

ID: PLP-SZY-2L

Cell Culture Media Bottle



- Square shoulder
- PET
- 500ml capacity
- Bulk packed
- Sterile



ID: PLP-PET-500-T-J

\$120.*
.00

BenchTop Essentials SALE Terms and Conditions



BenchTop Essentials SALE Specials are subject to change and only available from 01/03/2026 to 30/05/2026. Not valid with any other special offer or discounts. Special prices are as quoted. All prices are in AUD and GST is applicable on all purchases. Local delivery fees apply. Further product information and specifications including procurement Terms and Conditions of Sale can be found at: www.pacificlab.com.au or by contacting sales@pacificlab.com.au

Images displayed on this publication are a graphic representation of the product only and may vary from the actual item itself.

Pacific Laboratory Products

-  /PacificLabP
-  @PacificLabP
-  Pacific Laboratory Products
-  Pacificlaboratoryproducts

Blackburn, VIC 3130
Free Call: 1800 723 405
Ph: (03) 9845 0300

e-mail: sales@pacificlab.com.au