




PLP Bio-Science Plasticware

ELISA Plates 8 Strip Well, Flat Bottom



- Made from high grade USP Class VI virgin polystyrene providing high protein binding of more than 600ng/cm².
- Excellent optical clarity with CV of less than 5%.
- Flat Bottom, 96 well non treated plates with low refractive index and high light transmittance.
- Detachable, smart and economic.

Part #:	PLP-ESP-96-D
Type:	8-strip Well, Flat Bottom
Material:	Hi grade PS
Colour:	Transparent
Packaging:	10 pcs / bag, 200 pcs / box
Case Size:	47.5 x 29.5 x 33 cm

