

PCR TUBE / PCR 8-STRIP TUBE

White



White PCR 8-strip tubes improve the sensitivity and consistency of qPCR data by preventing fluorescence from refracting out of the tube

Polypropylene, Thin and Uniform Tube Wall, Heated Evenly



Advanced processing technology, flat cover with very low light loss



Good sealing to avoid cross-contamination



Transparent



Product Name	Cat. No.	Spec.
PCR 8-strip Tube, 100 μ L, with Cap	PCR-0801-FC	125 sets/box, 10 boxes/ctn
PCR 8-strip Tube, 200 μ L, with Cap	PCR-0802-FC	125 sets/box, 10 boxes/ctn
White PCR 8-strip Tube, 100 μ L, with Cap	PCR-0801W-FC	125 sets/box, 10 boxes/ctn
White PCR 8-strip Tube, 200 μ L, with Cap	PCR-0802W-FC	125 sets/box, 10 boxes/ctn
PCR Tube (Flat Cap)	PCR-02-FC	1000 pcs/box, 10 boxes/ctn