

# Labnet Pipette Selection Guide

Labnet International, Inc.

Labnet International offers a wide selection of manual and electronic pipetting products to choose from. Select your liquid handling product features below.



Products	Volume Settings	Channels	Volume Range	Pipette Types	Electronic Operation	Finger hook	Autoclavable	Multi-Dispensing	In-Lab Calibration
BioPette™ Plus	Push button Adjustment knob	Single	0.1-10,000 µl	n/a	✗	✓	✓	✗	✓
BioPette™ Plus Multichannel Pipettes	Push button Adjustment knob	8 or 12	1-300 µl	n/a	✗	✓	✓	✗	✓
BioPette™ Single Channel Pipettes	Push button Adjustment knob	Single	0.1-10,000 µl	n/a	✗	✗	✓	✗	✓
BioPette™ Multichannel Pipettes	Push button Adjustment knob	8 or 12	1-300 µl	n/a	✗	✗	✓	✗	✓
Excel™ Electronic Pipettes	Programmable	Single	0.1-1.0 µl	n/a	✓	✓	✓	✓	✗
Excel™ Multichannel Electronic Pipettes	Programmable	8 or 12	0.5-1,200 µl	n/a	✓	✓	✓	✓	✗
Labpette™ FX Fixed Volume Pipettes	Fixed volume	Single	2-1,000 µl	n/a	✗	✓	✓	✗	✓
Labpette™ R Repeating Pipettes	Click stop	Single	1-5,000 µl	n/a	✗	✗	✗	✓	n/a
FastPette™ V2 Pipette Controllers	Finger switches	Single	0.5-100 ml	Glass or plastic	✓	n/a	✓	✓	n/a
FastPette™ Pro Pipette Controllers	Finger switches	Single	0.5-100 ml	Glass or plastic	✓	n/a	✓	✓	n/a

Manual



Controllers



Repeating



Electronic

