

96 & 384 Well Microplates

Compatible with Roche LightCycler® 480

Real-Time PCR System




The new Axygen brand microplates for the Roche LightCycler® 480 Real-Time PCR system are a perfect fit to maximize the potential of the LightCycler 480. The consistent thin wall design presents optimized assay conditions, while plates are manufactured from white polypropylene to magnify signal intensity and prevent “cross-talk” between wells. Printed black alphanumeric make it easier to locate wells and orient the plate. All plates are free of any detectable Rnase/Dnase and endotoxins (pyrogens).

- ▶ Optimized for the plate-based Roche LightCycler 480 Real-Time PCR System
- ▶ Available with optional bar coding and sealing film
- ▶ Centrifugation to 4,000 RCF (Relative Centrifugal Force)
- ▶ Cut-away notches for perfect plate positioning
- ▶ Consistently thin walls for reduced cycle times and reproducible results
- ▶ Automation friendly SBS footprint 384 well plate


96 Well Microplate

- ▶ For reaction volumes of 10-100µL
- ▶ Half skirt design
- ▶ Also available in Clear



384 Well Microplate

- ▶ For reaction volumes of 5-20µL
- ▶ Full skirt design
- ▶ Also available in Clear



Ordering Information

Cat.No.	Description	Qty/Unit	Qty/CS
PCR-9 6-LC480-W	96 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White (with UC-500 Sealing Film)	10 Plates & Films	5 Units
PCR-96-LC480-W-NF	96 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White	10 Plates	5 Units
PCR-96-LC480-W-BC	96 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White, Bar Coded(with UC-500 Sealing Film)	10 Plates & Films	5 Units
PCR-96-LC480-C-NF	96 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, Clear	10 Plates	5 Units
PCR-384-LC480-W	384 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White (with UC-500 Sealing Film)	10 Plates & Films	5 Units
PCR-384-LC480-W-NF	384 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White	10 Plates	5 Units
PCR-384-LC480-W-BC	384 Well Microplate Compatible with Roche LightCycler 480 Real-Time PCR System, White, Bar Coded (with UC-500 Sealing Film)	10 Plates & Films	5 Units
UC-500	UltraClear Pressure Sensitive Plate Sealing Film	100 Films	5 Units

CORNING | FALCON® AXYGEN® PYREX®

Distributed by PACIFIC LABORATORY PRODUCTS

Australia Head Office

PO BOX 465
Blackburn, VIC 3130

Free Call: 1800 723 405
Ph: (03) 9845 0300
Fax: + 61 3 9845 0350

e-mail: sales@pacificlub.com.au

www.pacificlub.com.au

/PacificLab @PacificLabP