

# Labnet Mini Microcentrifuges

- Quick release rotor system
- Electronic braking
- Dual lid switch
- Compact, rugged design



## SPECIFICATIONS

Maximum Speed	6,000 rpm / 2,000 x g
Maximum Capacity	8 x 1.5/2.0 mL or 4 x 0.2mL strips
Dimensions	6.8 x 5.8 x 4.8 in / 17.2 x 14.7 x 12.2 cm
Electrical	AC 100 - 240V, 50/60Hz

## ORDERING INFORMATION

C1601-B	Labnet Mini Microcentrifuge, blue lid, 100-240V	
C1601-P	Labnet Mini Microcentrifuge, purple lid, 100-240V	
C1601-R	Labnet Mini Microcentrifuge, red lid, 100-240V	
C1601-G	Labnet Mini Microcentrifuge, gray lid, 100-240V	
6770-RT	Rotor for 8 x 1.5/2.0mL tubes	
6770-RTS	Rotor for 4 x 8-strip tubes	
6770-RTA.20	Individual adapters for 0.2mL tubes, fits in rotor 6770-RT, pack of 8	
6770-RTA.25	Individual adapters for 0.25mL tubes, fits in rotor 6770-RT, pack of 8	
6770-RTA.5	Individual adapters for 0.5mL tubes, fits in rotor 6770-RT, pack of 8	