

# BenchTop Centrifuge

## Compact Table Top Centrifuge from Pacific Laboratory Products

**12**  
Months  
Warranty

Power  
Ready  
Australian  
standards

*Solid, easy to use and reliable centrifuge  
ideal for Clinical Chemistry, Specimen  
Reception and Blood Collection Centres*



- Intuitive microprocessor control panel
- Stainless steel, swing out rotor
- Brushless induction drive motor
- Safety lid-lock mechanism
- Stainless steel reinforced housing
- Automatic lid opening system
- Range of tube adaptors available


**PATHO1**




- Stainless Steel Housing
- Swing Out Rotor

### Specifications


Max Speed	4000 RPM
Max RCF	2220 x g
Swing Rotor	8 x 15 ml
Time Range	0 ~ 99 min
Power Supply	AC 220V-240V, 50Hz, 5A
Noise	≤ 65 db(A)
Weight	15kg
Dimensions	330 x 430 x 290mm (W x L x H)



PATHO1AD1



PATHO1AD2



PATHO1AD3