



## IMPORTANT

### Safety Requirements:

- \* After an installation, the purchaser is responsible to register the autoclave (pressure vessel) with the local Occupational Health and Safety authority (Work Safe).
- \* Autoclaves are considered pressure vessels and so require to be certified that they are safe and meet Australian standards. This task is done by a third party certification company in Australia (unrelated to Pacific Laboratory Products).
- \* Autoclaves also require scheduled maintenance and must be performed following strict manufacturer's recommendations.
- \* Be advised that annual inspections are undertaken by external (unrelated to Pacific Laboratory Products) authorized consultants and may incur fees.

If you require registration documents to support the design of the autoclave or would like to know more about these units, simply contact one of our sales representatives at: [sales@pacificlab.com.au](mailto:sales@pacificlab.com.au)

Regards,

The team at Pacific Laboratory Products.