
 Learn more! 
 Watch the Video

Hit the play button!



ZymoSnap ALP | Alkaline Phosphatase Test



ZymoSnap ALP rapidly measures alkaline phosphatase activity in dairy products to verify pasteurization. The test has the capability to detect below the regulatory limit (350mU/L) down to 100mU/L. The device features an all-in-one design that provides simplicity, convenience, and repeatability. In addition, the clear definitive results avoid delays and additional costs of retesting as seen with other ALP test methods. ZymoSnap ALP has been validated for a variety of milks, flavored milks, creams, and more.

Benefits

- All-in-one device
- 15 second read time
- Large sample throughput (30 tests per hour)
- 10 minute assay
- Suitable for milk products such as flavored milks, cream, yogurts, shakes
- Compatible with ENSURE luminometer
- 12 month shelf life at refrigerated temperatures (2-8°C) *
 (* From date of manufacturing).


Sensitivity

- Detects down to 100mU/L
- Clear, reliable results even at low levels



Positive Control Kit for ZymoSnap ALP Alkaline Phosphatase Test

ZymoSnap ALP Procedure



Part Number	Description		Quantity
HG-ZS-ALP-100	ZymoSnap ALP: Alkaline Phosphatase Detection Test	Shop	100
HG-ZSALP-PC	Zymosnap Alkaline Phosphatase Positive Control	Shop	100