

Nylon Mesh Cell Strainers



Made with Polypropylene frame

- ✓ Designed for Tissue Culture Labs
- ✓ Ideal for filtration of various Cells
- ✓ Strong Nylon Mesh for better performance
- ✓ Compatible with most brands of 50ml tubes
- ✓ Gamma radiation sterilised

40µm



ID: 15-1040

70µm




ID: 15-1070

100µm



ID: 15-1100



Dimensions (mm)

A = O.D. 35.7	E = 14.4
B = O.D. 25.7	F = 1.5
C = O.D. 20.6	G = 2.8
D = O.D. 7.8	H = O.D. 25.7



• Color coded for easy identification

• Non-pyrogenic

• Individual wrapped

• Moulded tab for easy handling

Code No	Product Information	Pack Size
15-1040	40µm Cell Strainer, Sterile, Individually wrapped, Purple	100 Pieces / Bag
15-1070	70µm Cell Strainer, Sterile, Individually wrapped, Orange	100 Pieces / Bag
15-1100	100µm Cell Strainer, Sterile, Individually wrapped, Green	100 Pieces / Bag

Distributed by



www.pacificlab.com.au

Australian Head Office

PO BOX 465
Blackburn, VIC 3130

Connect with us:

Free Call: 1800 723 405

Ph: (03) 9845 0300

Fax: + 61 3 9845 0350

e-mail: sales@pacificlab.com.au

 /PacificLabP  @PacificLabP