


## 50µl 72 Rotor Disc Qiagen, Rotor-Gene Compatible

- For use with Corbett & Qiagen Rotor-Gene instruments
- Contamination free – Free from human DNA, DNase, RNase and PCR inhibitors
- Virgin Resin – Only virgin polypropylene used
- Optically clear tubes
- Background fluorescence & evaporation tested
- Requires heat film



Part Number	Description	Pack Size
981303-24-4-SSP	50µl, 72 Rotor-Gene Rotor-Disc, Polypropylene (PP)	Box of 24

## Sealing Film for Rotor Disc

- To be used with 72 and 100 Rotor Discs
- Heat sealable, does require heating device
- Packed in tamper-evident re-sealable bags



Part Number	Description	Pack Size
981601-1-0-SSP	Rotor-Disc, Heat Sealing Film, For 100 and 72 Rotor-Discs	Bag of 60